

WILD M7A M7S

WILD
HEERBRUGG

Instructions for use
Mode d'emploi
Bedienungsanleitung
Instrucciones para el uso

Instructions for use

I. Assembly

1. Stands
 - 1.1 Swinging-arm and table-clamp stands
2. Optics carrier
 - 2.1 Stereoscopic and vertical-beam observation
3. Binocular tubes
4. Objectives
5. Eyepieces

II. Illuminators

6. Incident light
 - 6.1 Low-voltage lamp 6 V/20 W
 - 6.2 Low-voltage lamp 6 V/10 W
 - 6.3 Ring illuminator
 - 6.4 Fluorescent surface illuminator
 - 6.5 Coaxial incident illuminator
 - 6.6 Vertical incident illuminator
7. Transmitted light
 - 7.1 Transmitted-light stage
 - 7.2 Transmitted light stand, bright field, circular
 - 7.3 Transmitted-light stand EB
 - 7.4 Transmitted-light stand, bright/dark field
8. Combined illuminators
9. Regulating transformer

III. Accessories

10. Stages
 - 10.1 Gliding stage
 - 10.2 Cup stage
 - 10.3 Stage carrier with mechanical stage
 - 10.4 Hand-focused stage carrier
11. Handrests
12. Filterholder
13. Measuring
14. Attachment for vertical and oblique observation
15. Polarisation
 - 15.1 Polariser on glass stage plate
 - 15.2 Rotatable Pol. stage
16. Tubes
 - 16.1 Double-iris diaphragm
 - 16.2 Discussion tube
 - 16.3 Drawing tube
 - 16.4 180° dual stereo attachment
 - 16.5 Phototube
 - 16.6 Beam splitter

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
IV. Operation

17. Interpupillary distance
18. Exit pupil position
19. Dioptic correction
20. Focusing
21. Zoom magnification changer
22. Exposure time calculator
23. Choosing the optics combination
 - 23.1 Working distance
 - 23.2 Total magnification
 - 23.3 Field diameter
 - 23.4 Resolution
 - 23.5 Depth of field
24. Care of the instrument

V. Appendix

Dimensions of the M7A Stereomicroscope

Table: Optical data

Do you know on what principles your WILD M7A or M7S Zoom Stereomicroscope is designed?

The stereomicroscope produces a spatial image of three-dimensional objects. Two separate beam paths which traverse a common main objective cause each eye to observe the object from a different direction. The two images are fused into a single plastic image.

Wild stereomicroscopes are designed with two separate optical trains, the axes of which run parallel to one another after passing through a common main objective. Because the intermediate image planes are parallel to the object plane, the whole field is uniformly sharp and the natural link between accommodation and convergence is retained.

The WILD M7S Stereomicroscope enables subjects to be viewed either stereoscopically or vertically from above, because the objective can be displaced to lie in the axis of the left beam path (figs. 1, 2).

This results in an improved correction of chromatic aberration in photography, in the parallax-free measurement of the horizontal distance between two points at different heights, and in the correct rendering of polarisation colours.

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
**Here's how to use your
WILD M7A or M7S
Zoom Stereomicroscope
properly:**

1. Slide your stereomicroscope down the column and secure it using the sleeve, the safety ring and the wing screw. Sections 1-2
2. Fit the binocular tube and additional tubes if required. Section 3
3. Insert the eyepieces. Secure them with the clamping screws. Section 16
4. Choose an additional objective if required. Fit it to the main objective using the screw. Set the working distance. Section 5
5. Mount the illuminator required and switch it on. Section 4
6. Look down the binocular tube. Displace the eyepiece tubes to give the correct interpupillary distance. Section 6-9
7. Correct your diopter setting by turning the eyelens mounts. Section 17
8. Position the object to be investigated under the objective. Use the focusing drive to bring the image into focus. Section 19
9. Select a low magnification in order to obtain an overall view of the object. Rotate the zoom magnification changer to the desired magnification. Section 20
10. Section 21

Your WILD M7A or M7S Zoom Stereomicroscope is now parfocally adjusted. The image remains sharp from the lowest magnification to the highest. You can inspect the object, right from the overall view down to fine details, without re-focusing. Make sure, though, that you have gone through the procedure for dioptric correction described in section 19.

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
Do you know the advantages of your WILD M7A or M7S Zoom Stereomicroscope?

The WILD M7A and M7S Zoom Stereomicroscopes are characterised by excellent flatness of field. Further, the multiple-coated optics result in first-class image transfer characteristics and consequently in a crisp image.

The human-engineered design ensures a comfortable and fatigue-free working position with conveniently-situated controls. A rich choice of modular accessories solves many problems.

The world-renowned Swiss quality of the Wild instruments guarantees the long-term untroubled functioning of your stereomicroscope and is supported by the technical advice and the servicing facilities available in the Wild agencies which are established in over one hundred countries.

I. Assembly

The WILD M7A and M7S Zoom Stereomicroscopes are of the same basic design and are equipped with the same superlative Wild optics.

The complete outfit comprises an optics carrier with drive housing, inclined binocular tube, main objective, eyepieces and optional additional objectives. These components are complemented by the appropriate stand, a suitable illuminator and the correctly-chosen modular accessories.

Both instruments produce an erect, laterally-correct, stereoscopic image which permits spatial subjects to be studied and manipulated with the convenience of a long working distance.

The WILD M7S Zoom Stereomicroscope possesses in addition a unique feature - the main objective can be displaced to lie coaxially in the left beam path, and consequently the subject is viewed vertically from above instead of at an angle.

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



4

1. Stands

Numerous incident- and transmitted-light stands are available for the WILD M7A and M7S Stereomicroscopes.

Incident- and transmitted-light stands

The normal stands for incident light are circular or rectangular.

The incident- and transmitted-light stands are supplied with a 350 mm long column. An inclinable 215 mm long column is available for the rectangular incident-light base (fig. 3). If the base and the column are found to have been packed separately, connect them with the two Allen screws, using the Allen key provided (fig. 1).

- Slide the safety ring down the column.
- Place the stage clips and the stage plate in position.
- Fit the optics carrier as described in section 2.

The stereomicroscope can be provided with a metal hood which fits on the circular incident- and transmitted-light stands, protects the instrument, and enables it to be easily carried.

Large incident-light stand

To inspect masks and wafers, large incident-light stands are available (fig. 2). Wafer slides and turrets, and mask holders, can be supplied to order.

- Bolt the stage carrier (29) to the baseplate from beneath.
- Slide the sleeve and the safety ring down the column.

This stand can also be equipped with a fibre ring lamp for observations in transmitted light; a carrier for it is to be placed in the baseplate.

Universal stand

The universal stand consists of a baseplate and a 450 mm- or 800 mm-long column which is to be bolted to the centre of the baseplate (fig. 4).

- Slide the sleeve and the safety ring (36) (flat side uppermost) down the column.
- Lower the adapter (37) so that the carrier rod (38) points upwards.
- Remove the safety screw (28).
- If appropriate, slide an adapter for column diameter 25 mm on to the carrier rod so that a lamp can later be added. Slide the stereomicroscope down the column.

The baseplate serves as a working surface and a transmitted-light stand can be placed upon it.

Large stands

Floor stands, wall stands and ceiling mounts are available for special applications; they have a cantilever of 140 cm.

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
1

1.1 Swinging-arm and table-clamp stands


Two models of swinging-arm stand and table-clamp stand are available for the WILD M7A and M7S Stereomicroscopes (figs. 1, 2). Each consists of a stable base or a clamp, a post (20), a horizontal arm (21) and an inclinable carrier rod (22).

The post and the horizontal arm of the large swinging-arm stand and of the large table-clamp stand are larger, so the radius of activity is extended (fig. 2).

The clamp of the table-clamp stands can be secured to table tops 20 mm–50 mm thick.

- Secure the cross-member (29) to the post with the clamping screw (24) and tighten the safety ring (23).
- Slide the horizontal arm into the cross-member and secure it with screw (25).
- Slide the carrier rod over the horizontal arm and secure it with the clamping screw (26).
- The clamping lever (27) enables the inclinable carrier rod to be fixed at any angle. The two clamping levers can be pulled out along their axes and reset at another position.
- Remove the end-stop (28).
- Mount the stereomicroscope and, where appropriate, the low-voltage lamp (see section 6) on the carrier rod. Replace the end-stop.

2

The following movements are possible with the swinging-arm and table-clamp stands:

- The horizontal arm can be swung out sideways.
- The carrier rod can be tilted
- The working distance is variable
- The horizontal arm can be extended.

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



2. Optics carrier

The M7A and M7S optics carriers have an integrated drive housing with maintenance-free ball-bearing slides and a 25 mm diameter bore for attachment to a stand.

- Slide the optics carrier down the column of the stand, to rest on the safety ring. The safety ring serves to protect the stereomicroscope, and also as a pivot for lateral movement.

Screw (4) enables the binocular tube, or an accessory tube, to be attached to the optics carrier.

M7A optics carrier

Position the M7A optics carrier over the stage and secure it with the clamping screw (6) (fig. 1).

When screw (96) is loosened, the stereomicroscope can be fully rotated and fixed in any position.

M7S optics carrier

To centre relative to the incident-light stand (fig. 2), swing the main objective into the axis of the optics carrier (see section 2.1).

- Tighten the securing screw (98).
- Move the changeover lever (97) to the left so that the instrument is against the end-stop.
- Loosen screw (98).
- Position the instrument over the centre of the stage plate.
- Tighten screw (98); the clamping screw (6) remains loose so that the instrument can be swung into the vertical-beam position.

To centre relative to the transmitted-light stand, it is recommended that the objective be swung into the left-hand beam path.

- Remove the left eyepiece and place a frosted glass over the empty tube.
- Displace the stereomicroscope laterally and adjust the light spot at low magnification so that the tube outlet is well illuminated.

2.1 Stereoscopic and vertical-beam observation

For stereoscopic observations, the main objective is positioned in the axis of the optics carrier (fig. 3) by activating the changeover lever (97) and simultaneously displacing the entire optics carrier with binocular tube as far as the end-stop. The optical beam path now coincides with the centre of the stage plate (fig. 4).

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
3. Binocular tubes

The standard outfit includes the 45° inclined binocular tube (fig. 1, right).

For outfits with transmitted-light stands or additional tubes, the inclined binocular tube of low viewing height is recommended (fig. 1, left).

The straight binocular tube (fig. 2, right) serves to ensure a comfortable working position when the instrument is inclined.


The low inclined binocular tube and the straight tube are also available as version T with built-in double-iris diaphragm and with drive for adjusting interpupillary distance.

A tube wedge enables the inclination of the binocular tubes to be increased or decreased by 15°.

The binocular tube with variable angle of observation provides for fatigue-free body positions with all types of instrument assembly and with any combination of accessories (fig. 2, left). The angle of observation can be steplessly varied by up to 90° in either direction from the vertical and can also be set for either of two viewing heights. The image automatically remains erect and laterally-correct.

- The binocular carrier tube can be removed from the optics carrier by loosening the clamping screw and can be fitted in either of two opposing positions.
- Fit it as required and wriggle it into position so that the locating screw engages the slot.
- Tighten the clamping screw.

Accessory tubes, such as the drawing tube (camera lucida), double-iris diaphragm or phototube, are fitted between the optics carrier and the binocular tube in the same way (see section 16).

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
4. Objectives

The optics carrier includes a built-in main objective of magnification 1.0 \times and working distance 91 mm (see table p. 71).

Working distance, field diameter and total magnification can be altered by securing additional objectives with factors 0.3 \times , 0.5 \times , 1.5 \times or 2.0 \times in front of the main objective (fig. 1).


The 0.3 \times additional objective, with its working distance of 265 mm, is mostly used with the universal stand or with the swinging-arm stand, or if the stereomicroscope is to be mounted on a machine.

The 0.5 \times , 1.5 \times and 2.0 \times additional objectives can be used on all incident- and transmitted-light stands.

Depending on the working distance of the additional objective, the instrument must be raised or lowered on the column of the stand. Because of the short working distance of the 2.0 \times objective, the safety ring must be removed (fig. 2). By loosening the screw (18), the outer mount of the 2.0 \times objective can be removed to provide more space for working on the object.

If the standard objective is used alternately with and without the 2.0 \times additional objective, the corresponding working distances can be set without changing the height of the stereomicroscope on the column:

- Secure the 2.0 \times additional objective.
- Wind the focusing drive completely down.
- Slide the instrument up the column until the object is visible in the field of view.
- Secure the stereomicroscope with the screw (6) or (98).
- Use the focusing drive to focus on the object.
- Remove the 2.0 \times additional objective.
- Use the focusing drive to set the optics carrier to the working distance of 91 mm, and focus on the object.

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



6. Incident light

6.1 Low-voltage lamp 6 V/20 W

The 6 V/20 W incident lamp is a powerful source of intense and very uniform illumination with a colour temperature of 3200 K.

The articulated lampholder is universally inclinable and can be secured in any position using the single centrally-positioned clamping knob (55). Support the lamp housing when loosening the clamping knob.

The lampholder can be attached directly (fig. 2), or with an adapter, to all posts of 25 mm diameter.

Fitting the bulb:

- Remove the rear part of the lamp housing.
 - Hold the 6 V/20 W halogen bulb in a cloth and push it into position. Do not touch the quartz envelope with the fingers.
 - Close the housing with the wire clip.
- Connect the lamp with the 7 V regulating transformer (see section 9) and switch it on.

A new bulb must always be recented:

- Pull out the condenser lens (56).
- Cause the light spot to fall vertically on a white paper about 25 mm away.
- Using screw (57), focus the image of the filament and its mirror image on the paper (fig. 4a).
- Displace the filament images with the horizontal movement screw (58) so that they are opposite (fig. 4b).
- Using the vertical movement screw (59), bring the two images into contact (fig. 4c).

When the condenser lens is fully pulled out, a sharply-defined illuminated field of diameter 16 mm is formed about 75 mm in front of it. By retracting the lens, the diameter can be increased to about 40 mm without changing the position of the lamp housing. A frequently-used mean position is denoted by a click stop.

Larger illuminated fields can be produced by displacing the lamp housing.

A heat-absorbing filter is built into the lamp. A diffusing filter in swing-out mount is attachable to one of the rods (60). Additional filter mounts are available, so that a total of four 32 mm diameter filters can be attached.

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



6.2 Low-voltage lamp 6 V/10 W

The 6 V/10 W low-voltage lamp is for illuminating opaque objects obliquely from above.

The 6 V/10 W low-voltage lamp can be provided with an inclinable lampholder (45) or with a clampable one (46).

The clampable lampholder can be secured in any position by means of a central clamping knob (47). When loosening this clamping knob, hold the lamp housing. Because the connectors at the two ends of the arm are identical, the arm can be attached to adapters and to the rotatable coupling (50) of the lamp. One or two low-voltage lamps can be attached to the mount of the main objective of the M7A Stereomicroscope (fig. 3):

- Unscrew and remove the ring (19).
- Slide either an adapter (48) and a spacing ring, or two adapters for double illumination, over the mount of the objective.
- Replace the ring (19).
- Connect the adapter, the lampholder and the lamp housing.
- Connect one or two low-voltage lamps to the 7 V regulating transformer (see section 9) and align relative to the object.
- The light spot can be individually focused by turning the lamp mount (51).

The lampholder can be attached to the column of a stand by means of an adapter diam. 25 mm, which will also connect it to a separate cast base (fig. 2), to the carrier rod of the swinging-arm stand, or to one of the two short posts on the rectangular stands for incident and transmitted light. This has the advantage that the relative positions of the illumination and the objective remain unchanged during focusing.

With the M7S Stereomicroscope, this arrangement prevents disturbance of the illumination when the main objective is swung out:

- Depending on the working distance, secure one or two adapters (49) beneath the safety ring or above the drive housing on the column of the stand.

Changing the bulb:

Note: Always use a cloth to hold the halogen bulb.

An arrow on the lamp mount and a line on the lamp housing indicate the position at which the two parts can be separated by pulling. The lamp mount clicks automatically into position again when refitted.

Two 50 mm diameter filterholders can be attached to each of the two pegs on the lamp housing and can be swung in as required.

A heat-absorbing filter is supplied with this lamp.

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


6.3 Ring illuminator

The ring illuminator, which produces light of daylight quality, is used for the shadow-free illumination of spatial objects. The intensity suffices for magnifications of up to 40 \times .

- Remove the securing ring (19).
- Slide the holder ring over the objective mount and secure it with the securing ring.
- Using the screw (42), secure the lamphousing at the desired height.
- Set the choke to the available voltage as follows:

Unscrew the fuse holder (41).

Pull out the cover plate of the voltage selector and replace it so that the appropriate voltage is visible in the cutout:

110 V position for 100–120 V supply

220 V position for 200–240 V supply

Tighten the fuse holder.

- Connect the power cable.
- Switch on the illuminator with the toggle switch (40).

Fuses:

for 110 V: 200 mA, 5 \times 20 mm (370 879)

for 220 V: 500 mA, 5 \times 20 mm (167 650)

6.4 Fluorescent surface illuminator

The fluorescent surface illuminator produces light of daylight quality and serves to light up large surfaces without producing hot-spots or shadows (fig. 3).

- Connect the illuminator to the power supply, taking account of the voltage available. Two models are available, for 110 V/60 Hz and for 220 V/50 Hz.

- Each lamp can be switched on separately. For starting without flicker, press the white button for about two seconds. This will also increase the life of the fluorescent tubes.

Switch off by briefly pressing the red button.

- The best lighting arrangement is obtained by placing the ballast unit behind the stereomicroscope stand, bringing the two illuminator arms past the sides of the instrument.

- The angle of incidence of the light and the distance from the specimen are adjustable by means of joints. The joints are separately adjustable by means of screws, in order to provide the degree of braking required.

- The four white fluorescent tubes can be replaced by four UV fluorescent tubes (4 W, $\lambda = 300\text{--}400$ nm, with a maximum at 356 nm, e.g. Sylvania F4T5 BLB).

To change the tubes, first slacken the two screws on the reflector to remove the screen (39).

Avoid looking directly at switched-on UV tubes.

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



6.5 Coaxial incident illuminator

Coaxial illumination is available for flat specular subjects such as integrated circuits or metal sections. The light from a 6 V/20 W halogen lamp is directed along the two beam paths. Built-in polarising filters eliminate stray reflections.

The M7S Stereomicroscope requires a filter-holder with quarter-wave foil in front of the objective (fig. 3), so that complete extinction does not result. The degree of extinction is varied by rotating the quarter-wave foil.

The lamp housing and the stereomicroscope together form a compact unit which fits between the optics carrier and the phototube (see section 16.5).

- Remove the binocular tube (see section 3).
- Fit the coaxial illuminator to the optics carrier so that the lamp housing is on the right and pointing backwards at 45°. Manoeuvre it into position, ensuring that the locating screw engages the guide, and tighten it.
- Position the binocular tube on the coaxial illuminator and secure it.

For photography with the coaxial illuminator, a quartz plate must be built into the camera. Detailed instructions for the use of the coaxial illuminator are to be found in instruction booklet M2 232.

Note: The coaxial illuminator introduces a magnification factor of 1.5x.

6.6 Vertical incident illumination

A prism enables the concentrated light beam from the 6 V/10 W low-voltage lamp to be directed vertically downwards on to the object so that depressions or holes can be illuminated (fig. 4).

- Attach the 6 V/10 W low-voltage lamp to the mount of the main objective (see section 6.2).
- Fit the mount of the prism to the ring of the main objective or to the mount of the additional objective.
- Align the prism (62) towards the column of the stand and secure it.
- Position the low-voltage lamp horizontally so that the light beam impinges on the prism.
- Focus the light beam and illuminate the object as well as possible by displacing the prism.

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

7. Transmitted light

7.1 Transmitted-light stage

Place the transmitted-light stage in the stage cutout (fig. 2).

Fit a glass stage plate and the two stage clips. The light source can be daylight or the 6 V/10 W low-voltage lamp. The lampholder can be attached to the post (89). Cause the light spot to fall on the tiltable mirror.

Slight deliberate disadjustment of the mirror will produce inclined transmitted illumination, which is very advantageous for observing small translucent objects.

The positioning of the optics carrier on the transmitted-light stand is described in section 2.

7.2 Transmitted-light stand, bright field, circular

In the circular transmitted-light stand (fig. 3) for bright field, the light from the 6 V/10 W halogen bulb is directed through the specimen and into the objective by means of a deviating mirror, a collector lens, and a condenser lens.

– Assemble the outfit as described in sections 1–5.

– The cable is packed underneath the stand; connect it to the 7 V regulating transformer (see section 9).

– Move the lever (30) to its limit and extract the guide (34) completely. Hold the 6 V/10 W halogen bulb with a cloth, insert it, and replace the guide.

The light spot is focused by displacing the guide with the lever (30) so that each field of view can be optimally illuminated. Since the illumination is very bright, it is however recommended that the entire glass stage plate (diam. 80 mm) be illuminated.

In front of the collector lens is a diffusing filter. Another 32 mm diameter filter can be substituted for it, as follows:

– Pull out the lamp cable.

– Loosen the three screws (33) on the underside of the base and remove the cover plate.




– Remove the spring clip from the filter-holder (32).

– Insert the filter and replace the spring clip.

– Replace the cover plate.

If even more light is required for special applications, a 6 V/20 W halogen bulb can be inserted.

The metal hood can be attached to the pair of slits (31).

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
7.3 Transmitted-light stand EB, bright field

A complete illumination system including a 6 V/20 W halogen bulb, a regulating transformer and a thermal fuse, is built into the transmitted-light stand EB (fig. 1). The specially-designed collector system permits optimal use of the light. This stand is particularly useful for photomicrography in transmitted-light bright field.

Further details about the operation of this stand and in particular about the use of the interchangeable collector lenses and the swing-in frosted filter are given in the instruction booklet M2 295.

7.4 Transmitted-light stand, bright/dark field



The transmitted-light stand for bright and dark field (fig. 3) enables observations and photography to be carried out in either bright field or dark field.

In the dark field technique, no light passes directly through the specimen into the objective (fig. 2); light reflected at the interfaces of components with differing refractive indices travels into the objective, however, so that a bright image on a dark background is seen and the contrast of weakly-coloured or transparent specimens can be increased. The built-in 12 V/100 W halogen bulb produces high light intensity at a colour temperature ideal for photomicrography.

It is connected to the 12 V regulating transformer, which simultaneously controls the fan.

The changeover from bright to dark field and back is carried out with the lever at the rear of the base.

Further details on the use of this stand are to be found in the instruction booklet M2 280.

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

8. Combined illuminators

Numerous objects simultaneously display many of their structures only when they are illuminated by a combination of lighting techniques.

Combined inclined and coaxial incident light (fig. 1) enables objects possessing both flat specular areas and uneven parts, such as hybrid circuits and mechanical components, to be inspected.

A 6 V/10 W or 6 V/20 W low-voltage lamp and the 6 V/20 W coaxial illuminator can be connected simultaneously to a 7 V regulating transformer (see section 9) which accepts a maximum loading of 50 W.

The combination of incident and transmitted light (fig. 2) is recommended for partly-transparent objects such as fibres, insects, and textiles. The inclined incident light emphasises details and structures of the surface. Two low-voltage lamps can be mounted on the transmitted-light stand EB and connected to the built-in transformer.

For photomicrography, the low-voltage lamps should be connected to a separate MTr28 transformer, in order that the colour temperature of the halogen lamp remains constant. The same is valid for the combination of the low-voltage lamp with the transmitted-light stand for bright and dark field.

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




9. Regulating transformer

The 0-7 V/50 VA regulating transformer (fig. 1) will accept the 6 V/10 W and 6 V/20 W low-voltage lamps, the circular stand for transmitted light with 6 V/10 W halogen bulb, the attachment for vertical and oblique observation and the coaxial illuminator with 6 V/20 W halogen bulb. A maximum loading of 50 W is to be connected to the two sockets (66) (fig. 2).

The 0-12 V/100 VA regulating transformer (fig. 3) is used with the transmitted-light stand for bright and dark field with 12 V/100 W halogen bulb and with the coaxial illuminator if a 12 V/20 W halogen bulb is fitted. It has three sockets which accept a maximum combined loading of 100 W (fig. 4).

- Set the voltage selector (67) to the voltage available:

Position 115 for 100-120 V

Position 220 for 200-240 V

- Connect the power cable to socket (65) and to the power supply.

Setting the limiting voltage:

The click stop 6 V or 12 V prevents overloading of the bulb. Turn the regulating knob (63) clockwise until it engages.

Slacken the screw on the metal ring (64).

Turn the regulating knob back to zero.

Connect the lamp to the socket (66). Turn the regulating knob until the voltmeter shows 6 V or 12 V, as appropriate.

Tighten the screw on the metal ring.

Changing the fuse:

Pull out the power cable.

Remove the screw (68) from the underside.

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




III. Accessories

10. Stages

The use of special stages simplifies work with Wild stereomicroscopes; they can be used on all stands.

- Remove the stage plate and stage clips supplied as standard (fig. 1).
- Move the stereomicroscope correspondingly higher on the 350 mm column and secure it with the wing screw and the safety ring.

10.1 Gliding stage

The gliding stage (fig. 2) can be quickly moved in any direction and rotated.

- Place the stage in the cutout in the base of the incident- or transmitted-light stand and secure it with the lever (81).
- Fit a glass or black/white stage plate and the two stage clips.

10.2 Cup stage

The cup stage (fig. 3) is essential for examining cultures in petri dishes which are held in position by a movable holder. Insects and plants are secured to the rubber surface with needles.

- Place the supporting ring in the stage cutout.
- Add the cup stage, which can be tilted to any position.
- The freedom of movement of the cup stage can be further increased by mounting it on the gliding stage (fig. 4).

10.3 Stage carrier with mechanical stage

The stage carrier with mechanical stage (fig. 5) facilitates the systematic scanning of specimens.

- Place the stage carrier in the stage cutout and secure it with the lever (81).
- Add the appropriate stage plate (glass or black/white, as appropriate).
- Secure the mechanical stage to the stage carrier with the Allen screw (82). The mechanical stage can also be used on the rotating Pol. stage.

The object is secured with the arm and is moved in an x- or a y-direction (27 mm x 97 mm) with the knobs of the mechanical stage.

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



10.4 Hand-focused stage carrier

This stage carrier enables the object to be brought into focus without changing the height of the tube. This is advantageous with outfits for television microscopy and with the drawing tube.

The stage carriers can be mounted on the column of the incident-light stand for bright or dark field and can be equipped with a glass or black/white stage plate, or with a special stage (fig. 1).

The range of movement is 50 mm.

11. Handrests

The individually-adjustable handrests ensure precise and fatigue-free work on the object (fig. 2).

A rubber ring on the base prevents sliding on the table surface. The upper surface of the handrest is insulated with a layer of paint. It is recommended that extension pieces be added if the handrests are used in conjunction with an additional special stage on the base of the stand.

12. Filterholder

A 75 mm x 75 mm filterholder, which accepts colour correction filters, is available for Wild stereomicroscopes.

- Fit the filterholder over the objective mount and secure it with a screw (fig. 3).

If the WILD M7S Stereomicroscope is used in vertical-beam position with coaxial incident light, a rotatable filterholder for 50 mm x 50 mm quarter-wave foils is necessary (see section 6.5).

A filter slide housing is also available. This will fit between optics carrier and binocular tube (fig. 4). The filter slide accepts two filters which can be readily changed. Ordinary gelatine filters can be cut into strips to match the filter hole size of 23 mm x 45 mm.

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
1

13. Measuring

Because of their high optical quality, the WILD M7A and M7S Stereomicroscopes are particularly suitable for precise measurements.

The M7S is to be used preferably in the vertical-beam position (see section 2.1). If measurements are carried out with the right eye, e.g. with the measuring eyepiece of the Digital Length-Measuring Set, the binocular tube can be rotated 180° (see section 3). Depending on the outfit, request a shorter for the incident- or transmitted-light stands if necessary.

For measurement work, adjustable WILD eyepieces can be furnished with a graticule (reticle) in mount (see section 5).

The graticules available have various scales or grids, or a crosshair.

Angles can be measured with the 10 × wide-field goniometer eyepiece, or with the rotating stage used for polarisation work. Carriers are available for Leitz measuring stages (fig. 2).


Further details of procedures used for measuring are given in the instruction booklet M2 275.

The WILD MMS 235 Digital Length-Measuring Set (fig. 1) (with printer) streamlines and accelerates the action of measuring.

The determinations of the calibration values and of the measurement data are microprocessor-controlled. After a single calibration procedure for each magnification, the length of the feature is constantly visible on a digital display.

Calculation errors and check measurements are superfluous.

Detailed information is to be found in the instruction booklet M2 270.

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
14. Attachment for vertical and oblique observation

This attachment enables three-dimensional objects to be viewed and photographed either vertically or at an angle of 45° without it being necessary to turn or tilt them (fig. 1).

The attachment can be attached to the mount of the main objective and equipped with an illuminator.

Further details are given in the instruction booklet M2 265.

15. Polarisation

Double-refracting substances such as many crystals, petrological thin sections and fibres, can be more easily seen and identified in polarised light.



The M7S Stereomicroscope offers particular advantages in the determination of polarisation colours. The vertical-beam position enables the correct colour to be obtained, because the observation direction is exactly perpendicular to the plane of the section.

The transmitted-light stands accept either the polariser on glass stage plate or the rotating Pol. stage. All polarising sets require an analyser in rotating mount, which is fitted to the objective mount.

If photographs are taken in polarised light, a quartz plate must be built into the shutter-piece.

15.1 Polariser on glass stage plate

- Place the polariser on glass stage plate, polariser downwards, in the cutout of the transmitted-light stand (fig. 3).
- Rotate the analyser using the lever (83) so that the field of view without specimen is as dark as possible.
- Place the specimen in position and focus on it.

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
86 screw to secure rotatable stage

87 milled ring to rotate compensator




The rotatable Pol. stage is equipped with a swing-out polariser.

If required, a first-order red compensator (sensitive-tint plate) and a mechanical stage can be added to it. The compensator facilitates the identification of weakly-birefringent specimens.

- To fit the compensator, remove the screw (88) from the underside of the stage.
- Push the compensator in its metal mount fully into the stage.
- Retighten the screw.
- Fit the mechanical stage using the screw (82).
- Remove the glass stage plate and the stage clips from the transmitted-light stand.
- Fit the rotatable Pol. stage so that the two centring screws (84) are symmetrically positioned in relation to the column.
- Secure the stage with the lever (81).
- Switch on the light.
- Swing in the polariser with the lever (85) and pull out the compensator slide.
- Look into the binocular tube and move the lever (83) until the field of view without specimen is as dark as possible.
- Secure the specimen with the arm of the mechanical stage and focus on it.

With the help of a crosshair in the eyepiece (see section 5) and of the two centring screws, centre the stage as follows:

- The point about which the stage rotates is adjustable and must be brought into coincidence with the axis of the objective, i.e. with the centre of the crosshair.
- Spin the stage in order to establish the point of rotation in the field of view.
- Using the screws (84), displace the point of rotation in the directions ab and cd (fig. 3). Successively repeat this procedure, simultaneously rotating the stage, until the point of rotation coincides with the centre of the crosshair.
- Push in the compensator mount and rotate it until the field of view appears red. The rotatable Pol. stage is calibrated and can be used for angular measurements (see instruction booklet M2 275).

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

16. Tubes

The WILD M7A and M7S Stereomicroscopes have a removable binocular tube, so that various accessories such as phototubes, the drawing tube or the double-iris diaphragm can be fitted.

- Remove the binocular tube (see section 3).
- Mount the accessory tube.
- Replace the binocular tube on the accessory tube.

It is advantageous to use the low binocular tube with these combinations.

16.1 Double-iris diaphragm

With the double-iris diaphragm, which fits between the optics carrier and the binocular tube (fig. 1) or the phototube A, the depth of field can be matched to the characteristics of the object (see section 23.5). This is particularly advantageous in photomicrography in monocular assembly, because no double-iris diaphragm is incorporated in the phototube A.

Double-iris diaphragms are however an integral part of the binocular tubes T, of the trinocular phototubes HU and HV, and of the beam splitters.

The double-iris diaphragm enables a satisfactory practical balance to be reached between resolution and depth of field for any given magnification.

When the diaphragms are partly closed, the light cone entering the objective is reduced after passing through the magnification changer.

The depth of field increases and the image brightness is reduced.

The iris diaphragms are operated with the lever (71).

Position 1 is the smallest aperture; position "OPEN" is the largest.

16.2 Discussion tube

The WILD M7A and M7S Stereomicroscopes can be expanded into an excellent training instrument by combining them with a discussion tube (fig. 2). Two observers see the same image simultaneously, so that both instruction and discussion are possible.

The discussion tube is connected to the universal stand by means of a drive housing.

The optics carrier is fitted beneath the discussion tube, and the two binocular tubes are fitted above it.

Further details are given in the instruction booklet M2 130.

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

16.3 Drawing tube

In a drawing, characteristic features can be emphasised, less important detail can be omitted, and features in various planes of the object can be depicted. The drawing tube can be used in either incident or transmitted-light specimens in daylight and in a comfortable sitting position. The user sees the object and the drawing simultaneously. The scale of the drawing can be chosen at will with the assistance of the magnification changer and the additional objectives.

The mirror is protected from damage by a quartz coating.

- Fit the drawing tube between the optics carrier and the binocular tube so that the mirror points to the left or to the right.
- Move the lever (72) to the white spot so that the drawing surface is imaged in the binocular.
- Move the objective (74) of the drawing tube so that both the drawing surface and the pencil are seen in focus. The focusing range is from 100 mm to infinity.
- Place the object in position, switch on the illumination, and bring the image into focus.
- It is recommended that the drawing surface be brightened by means of a separate table lamp, e.g. with a Wild 25 W mains lamp on cast base.
- The red screw (73) releases the mirror which can be removed and turned in order to mirror in the images of documents or objects in various positions. Comparisons can thus easily be made.

If individually-made scales or grids are mirrored in and photographed with the object, measurements of number, length and area can be carried out.

The hand-focused stage carrier (see section 10.4) facilitates focusing, because the object can be brought into focus without changing the chosen image scale.

16.4 180° dual stereo attachment

The 180° dual stereo attachment enables two people to work simultaneously at a WILD M7A Stereomicroscope.

- Slacken the clamping screw (96) and rotate the optics carrier by 90°.
- Retighten the clamping screw.
- Secure the dual stereo attachment with the clamping screw (45) in the same manner as a binocular tube.
- Mount the binocular tubes for the two observers.

The dual stereo attachment has an incorporated double-iris diaphragm and a beam splitter (50%/50%) with two lateral outlets for photography, filming, TV or dual observation (see section 16.6).

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
16.5 Phototube

The WILD M7A and M7S Stereomicroscopes are ideal for photomicrography with one of the Wild Modular Photomicrographic Systems (fig. 2), for television and for cinematography.

Phototube A is for a monocular camera assembly; Phototubes HU and HV are for a trinocular camera assembly.

- Position the phototube over the left or the right beam path of the optics carrier and secure it.

The vertical-beam position is to be preferred for the M7S Stereomicroscope, under which circumstances the phototube with camera must be mounted over the left beam path.

- Accessories such as the coaxial incident illuminator, the double-iris diaphragm and the filter slide housing are to be fitted between the optics carrier and the phototube.
- Equip the HU and HV Phototubes with a binocular tube.

The transmitted-light stands EB and bright/dark field, and the 6 V/20 W incident lamp, are particularly appropriate for photomicrography on account of their high light output, homogeneous illumination and ideal colour temperature.

Camera for Wild Modular Photomicrographic System

To fit a Wild shutterpiece, first place an eyepiece tube (77) with a fixed eyepiece (9) in the phototube.

Detailed information is available in the instructions for your Wild Modular Photomicrographic System. -

Single-lens reflex camera


Commercially-available single-lens reflex cameras can be connected to Wild stereomicroscopes via a phototube with eyepiece tube (77), fixed eyepiece, connecting tube, 0.32x camera objective, adapter and T-mount.

TV and cine cameras

A TV camera or cine camera can also be very easily attached to a phototube.

For single-stage imaging, the connecting member is an adapter with C-mount (78).

For two-stage imaging, the connecting members consist of an eyepiece tube (77) with fixed eyepiece and an objective $f = 50\text{ mm}$ with C-mount.

WILD M7A Zoom Stereomicroscope with TV outfit
WILD M7S Zoom Stereomicroscope with WILD MPS 45/51 S Photoautomat

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
M7S con Photoautomat WILD MPS 45/51 S

16.6 Beam splitter

The beam splitter is the component to which the photo-, cine- and TV tubes and the tubes for a second observer are attached. Two versions of the beam splitter, with different beam partition characteristics, are available:

- for dual observation and television (50% to the binocular and 50% laterally).
- for photo and cine (20% to the binocular and 80% laterally).

The double-iris diaphragm, which is built in, influences the beams in the binocular and in both of the lateral outlets. The latter can be used simultaneously for different purposes (fig. 1).

- Fit the beam splitter between optics carrier and binocular tube so that the diaphragm lever (71) points towards the observer.
- The phototube, the cine/TV tube and the tubes for a second observer have a self-centring connection to the dovetail mount of the beam splitter.

Phototube

The phototube (fig. 2) can be used with commercially-available single-lens reflex camera. A T-mount (100) attaches these to the $f = 250$ mm or $f = 350$ mm adapters.

Cine/TV tube

The cine/TV tube (fig. 3) can be used for TV and cine cameras which accept a C-mount (101). The four possible engagement positions of the cine tube enable the camera to be positioned where it is least in the way.

A format-indicating graticule (reticle) for the phototube and for the cine/TV tube carries the format limits for 35 mm film for adapter tube lengths of 250 mm and 350 mm respectively. The graticule fits into the adjustable wide-field eyepieces and into those for spectacle wearers (see section 5).

Tube for second observer



This tube is available in a short (fig. 1) and in a long version.

Because the tube can be rotated about two axes, it can be positioned at any angle.

The knurled sleeve (104) enables the position of the image to be corrected.

Stereo attachment for second observer

The advantages of this accessory are the individual positioning in any direction and the comfortable binocular observation facility for the second observer.

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

IV. Operation

17. Interpupillary distance

The interpupillary distance can be individually set on all binocular tubes by displacing the two tubes.

- inclined binocular tube: 47 mm–76 mm.
- low inclined binocular tube, and straight binocular tube: 52 mm–76 mm
- binocular tube with variable angle of observation: 52 mm–76 mm; setting can be retained with screw (7).

If the interpupillary distance has been correctly set, a single circular field is seen, which results from the fusion of the two separate fields.

18. Exit pupil position

The correct distance of the eye from the eyepiece is important in the quality of the image. The further the eye is from the correct position, the greater the loss of the peripheral part of the field, i.e. the field diameter is reduced and the marginal areas become dark.

The correct eye-to-eyepiece distance is 12 mm for the adjustable wide-field eyepieces and 22 mm for the eyepieces for spectacle wearers. At the correct position, the exit pupil, all of the ray bundles emerging from the eyepiece are combined.


The exit pupil is located by moving the eye slowly towards the eyepiece until the whole field of view is seen sharply defined.

The 22 mm exit pupil of the eyepiece for spectacle wearers enables spectacles to be worn. The flexible eyecups (12) prevent scratching of the spectacles. They can be folded back, if the entire exit pupil distance is required.


In order that the eyepiece for spectacle wearers may also be used by non-spectacle wearers, spacing rings of 6 mm and 9 mm thickness respectively (17) are available, with which the observer can individually set the required degree of eye contact.

- Screw the required spacing rings into position between the eyecups and the eyepieces (fig. 5). The flexible eyecups ensure comfortable working.

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

19. Dioptric correction

The observer must correct any ametropia on both eyepieces. The correct diopter settings are essential if parfocality is to be maintained over the whole zoom range.

- a Set the stereomicroscope to the working distance of the objective (see section 23.1)
- b Set the interpupillary distance (see section 17)
- c Set the two eyelens mounts to zero
- d Open the double-iris diaphragm if used.

- 1 Turn to the highest zoom position (position 31).

Focus on a flat test specimen (e.g. a cross on paper) with the focusing drive.

- 2 Turn to the lowest zoom position (position 6)

Turn the eyelens mounts fully out in the '+' direction without looking into the eyepieces. Turn each eyelens mount individually in the '-' direction until the specimen appears sharp. While doing this, close the other eye.

The setting can be retained on eyepieces for spectacle wearers.

- 3 Set the magnification changer to the highest position.

Insert the sharpness of focus and correct it with the focusing drive if necessary. The specimen should now remain in focus at all zoom positions without any further refocusing.

Otherwise the procedure must be repeated. The diopter values set remain the same irrespective of the objective or eyepiece used.

If a graticule is used in the eyepiece, it must be in position while the above adjustments are being carried out, because it affects the diopter setting.

The procedure for dioptric correction with the graticule is described in the instructions for the photomicrographic systems and for measuring equipment.

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
20. Focusing

The bilateral focusing knobs (1) ensure comfortable operation. The superior intermediate optics of the WILD M7A and M7S Stereomicroscopes ensure that the image remains sharp over the whole zoom range, provided that the dioptric settings are correct and that the working distance is not changed.

Features nearer to, or further from, the principal plane of the object can be imaged by refocusing. The focusing range of 50 mm enables even relatively large objects to be brought into focus without raising the stereomicroscope on the column. The ease of movement of the focusing drive is adjustable to the loading by rotating the two drive knobs in opposite directions.


21. Zoom magnification changer

The advanced intermediate optics of the stereomicroscopes enable zoom positions to be set steplessly within the range 1:5 while the working distance remains constant. When the milled ring (2) is rotated, the distances between the components of the lens combinations are changed, so that the image remain sharp at all zoom positions.

The magnifications obtained with 10× eyepieces (from 6× to 31×) are inscribed on the milled ring. The black arrow on the optics carrier indicates the position (e.g. "11" in the illustration).

The screw (7) enables the magnification set to be retained.

The upper and lower limits of the zoom range are defined by stops.

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
22. Exposure time calculator

Beneath the zoom magnification changer are two scales which enable identical exposures to be obtained at different zoom positions.

The red scale (91) can be displaced, and enables an exposure factor to be determined.

Example:

- Engage zoom position "6".
- Meter the exposure time and take the photograph.
- Displace the red scale (91) so that the arrow with figure "1" is opposite the figure "6" on the black scale (92).
- Engage the required magnification position, e.g. "22".
- Look for the figure "22" on the scale (92).
- Read off the exposure factor "2" on the red scale (91) and expose correspondingly longer.

23. Choosing the optics combination

The choice of optics combination is influenced by the working distance required, the total magnification obtainable, and the appropriate field of view. Resolution and depth of field are also dependent on the optical components chosen.


The table (page 71) gives field diameters and total magnifications for zoom positions 6, 10, 20 and 31 in conjunction with the adjustable wide-field eyepieces 10 \times , 15 \times , 20 \times or 32 \times and the objective combinations 1.0 \times , 0.3 \times , 0.5 \times , 1.5 \times and 2.0 \times .

23.1 Working distance

The working distance can only be changed by substituting or adding an objective of different focal length. The higher the objective magnification, the shorter the focal length and the working distance.

Additional objective 0.3 \times	with working distance 265 mm
Additional objective 0.5 \times	with working distance 154 mm
Main objective 1.0 \times	with working distance 91 mm
Additional objective 1.5 \times	with working distance 38 mm
Additional objective 2.0 \times	with working distance 25 mm

To fit the additional objectives, see section 4.

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

23.2 Total magnification

The imaging within a stereomicroscope takes place in two stages.

In the first stage, the objective together with the magnification changer forms the intermediate image of the object.

In the second stage, the intermediate image is further magnified by the eyepiece.

The total magnification is the product of the magnifications of eyepiece and objective, and is calculated as follows:

$$\frac{M_{ob} \times M_{oc} \times p \times t}{10} = M_{total}$$

Example: $\frac{1.5 \times 15 \times 20 \times 1.5}{10} = 67.5 \times$

- M_{ob} magnification of objective
- M_{oc} magnification of eyepiece
- p zoom position
- t tube factor (e.g. 1.5× for coaxial illuminator)
- M_{total} total magnification

The magnification range can be extended by using additional objectives. The 0.3× and the 2.0× objectives produce lower and higher ranges respectively.

The 10×, 15×, 20× and 32× eyepieces also influence the magnification, so that total magnifications from 1.8× to 198× are attainable (see table p. 71).

23.3 Field diameter

The field of view is the circular area of the specimen which can be seen down the microscope (fig. 2). It is influenced by objective and eyepiece characteristics and decreases with increasing magnification, e.g. 117 mm at 1.8× and 2.1 mm at 124× (see table p. 71). The field diameter is calculated as follows:

$$\varnothing = \frac{\text{field number} \times 10}{M_{ob} \times p \times t}$$

The field numbers are inscribed on the eyepieces (e.g. 10×/21; field number: 21).

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



23.4 Resolution

The resolving power of a stereomicroscope is a measure of its ability to distinguish fine detail, i.e. to separate two closely-spaced points.

Resolution and contrast depend on the objective and on the zoom changer position. Resolution is proportional to the angle of aperture (fig. 1). Powerful objectives of short working distance thus give high resolution, and this increases with the zoom position.

Details not resolved by the objective cannot subsequently be rendered visible by using strong eyepieces.

A given total magnification is attainable either by using a weak objective coupled to a strong eyepiece, or vice versa, e.g. 62 \times with a 1.0 \times objective, 20 \times eyepieces and zoom position 31, or with a 2.0 \times objective, 10 \times eyepieces and zoom position 31.

The second combination is preferable if the working distance is unimportant.

An additional eyepiece magnification should be so chosen that all details resolved by the objective are recognisable. No advantage is gained by an excessive secondary magnification, but an inadequate secondary magnification means that the potential of the objective is not being fully exploited.

23.5 Depth of field

Depth of field depends on the objectives, the eyepieces and the tube factor (see our brochure "Practical Photomacro- and Micrography", M3 300e.)

Depth of field decreases with increasing magnification.

The user can use these laws to match the appropriate optics combination to the characteristics of the specimen.

For larger objects with coarse structures, the double-iris diaphragm can often be used to improve depth of field at the expense of resolution (see section 16.1).

Fig. 3: Double-iris diaphragm open
Fig. 4: Double-iris diaphragm partially closed

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
Fig. 4: Diafragma enclavado

24. Care of the instrument

Wild Stereomicroscopes are high-performance precision instruments of the highest optical and mechanical quality. The security offered is correspondingly high:

We guarantee the quality of each of our instruments. Our guarantee covers all faults in materials and manufacture. It does not, however, cover damage resulting from careless or improper handling.

It is evident that a valuable optical instrument merits considerate treatment. It will then fulfil its functions for decades with the same precision.

Dust and dirt have a very bad effect on the image quality.

- Use the dust cover or the metal hood, therefore, during pauses in work.
- Eyepieces, tubes without eyepieces, and the optics carrier when without tube, should always be protected against dust by means of the covers provided.
- Dust on optical components is best removed with a blow-bulb or with a soft brush.
- Tenacious marks on objectives and eyepieces can be removed with a soft lint-free cloth dipped in pure alcohol. Special cloths for cleaning optics are also available.
- Never dismantle optical systems unless instructions for doing so are given in the instruction booklet.
- Protect optical instruments from damp, fumes, corrosive substances, and acid or alkaline vapours. Do not store chemicals nearby.
- The mechanical parts of the instrument can be cleaned with a soft brush or with a soft cloth.
- The drive housing (with ball-bearing slide) is maintenance-free.
- Never grease, oil or dismantle mechanical parts unless instructions for doing so are given in the instruction booklet.

If the instrument no longer functions satisfactorily, call in a specialist or your local Wild agent, or contact the factory.

Modifications resulting from technical developments may be made in the interest of our customers. Thus, illustrations, specifications and catalogue references are not binding and are subject to change without notice.

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© Wild Heerbrugg Ltd
Wild Heerbrugg Ltd
Leica Mikrosysteme Vertrieb GmbH